

## MOULD INCUBATOR/ BIOCHEMICAL INCUBATOR



SPX-B-Z / MJX-B-Z  
Series

Shanghai Boxun Medical Biological Instrument Corp

Add: No. 599 Zhongqiang Road, Maogang Town,  
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com



www.sh-bx.com

## MOULD INCUBATOR/BIOCHEMICAL INCUBATOR »

### SPX-B-Z / MJX-B-Z Series

#### Product Features

- The shell is made of cold-rolling steel with static spray plastics, the working chamber is made of mirror stainless steel, the shelves can be adjusted optionally;
- PAMMA I control system;
- With timing functions;
- Lighting device in working chamber (mould incubator is equipped with UV light);
- Large viewing angle heat preservation vacuum toughened glass;
- Power supply recovery function, data will not be lost when power failure and system halted.



Model	SPX-100B-Z MJX-100B-Z	SPX-150B-Z MJX-160B-Z	SPX-250B-Z MJX-250B-Z
Volume [L]	100	150	250
Type	SPX Series Biochemical Incubator/MJX Series Mould Incubator		
Internal Dimensions W×D×H [mm]	450×380×585	500×335×955	505×450×1100
Housing Dimensions W×D×H [mm]	580×570×1070	595×580×1232	615×620×1597
Temp. Range [°C]	5-50	5-50	5-50
Temp. Resolution [°C]	0.1	0.1	0.1
Temp. Fluctuation [°C]	±1°C (Heating), ±2°C (Cooling)		
Temp. Uniformity at 37°C [°C]	±2		
Power [W]	300	300	360
Timer Range [h]	0-999h		
Number of Shelves (std./max.) [pc]	2/3	3/5	3/7
Load per Shelf [kg]	15		
Net Weight/Gross Weight [kg]	62/73	79/91	100/118
Package Size [mm]	680×670×1240	730×680×1400	720×720×1800