

VACUUM OVEN



BXK Series



BXK-S Series

Shanghai Boxun Medical Biological Instrument Corp
上海博迅医疗生物仪器股份有限公司

Add: No. 599 Zhongqiang Road, Maogang Town,
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com



www.sh-bx.com

Vacuum Oven >>

BXK Series

Product Feature >

- The shell is made of cold rolled steel plate, the surface is electrostatically sprayed, the inner is made of stainless steel;
- Microcomputer single-chip temperature control system, temperature control, timing, over temperature alarm functions
- Timing and timing functions ;
- Temperature can be controlled and displayed independently for each heating shelf;
- Integrated design, standard vacuum pump placed inside the noise reducing cabinet;
- The door can be adjusted, a silicone rubber door seal, toughened, bulletproof double glass door;
- Power supply recovery function, data will not be lost when power failure and system halted;
- Optional: Oil filter.

BXK-S Series

Product Feature >

- The shell is made of cold-rolling steel with static spray plastics, the working chamber is made of stainless steel with arc angle design which is convenient to clean;
- Microcomputer control system with temperature and time control and over temperature alarm, LCD screen;
- Timing and timing functions ;
- Temperature can be controlled and displayed independently for each heating shelf;
- Integrated design, standard vacuum pump placed inside the noise reducing cabinet;
- The door can be adjusted, a silicone rubber door seal, toughened, bulletproof double glass door;
- Power supply recovery function, data will not be lost when power failure and system halted;
- The vacuum degree can be displayed and controlled directly;
- Optional: Oil filter, intelligent touch screen controller.



Model	BXK-30(S)	BXK-50(S)	BXK-100(S)	BXK-250(S)	BXK-500(S)
Volume [L]	30	53	100	250	500
Temp. Range [°C]	RT+10~250	RT+10~250	RT+10~250	RT+10~250	RT+10~250
Temp. Fluctuation [°C]	±1	±1	±1	±1	±1
Vacuum Degree [Pa]	133	133	133	133	133
Power [W]	450	1200	1400	2800	4000
Power Supply [V, Hz]	220±10, 50/60	220±10, 50/60	220±10, 50/60	220±10, 50/60	380±10, 50
Timer Range [min]	0~9999	0~9999	0~9999	0~9999	0~9999
Internal Dimensions W×D×H [mm]	320×320×300	415×370×345	450×500×450	600×600×700	630×810×845
Housing Dimensions W×D×H [mm]	480×530×1015	575×550×1060	620×640×1165	770×780×1515	800×990×1660
Number of Shelves [pc]	1	2	2	3	4
Net Weight [kg]	84	109	129	228	442
Gross Weight [kg]	98	122	148	251	475
Package Size [mm]	610*580*1240	685*645*1275	720*775*1420	870*880*1650	1090*970*2060