

## VACUUM OVEN



BZF Series

Shanghai Boxun Medical Biological Instrument Corp  
上海博迅医疗生物仪器股份有限公司

Add: No. 599 Zhongqiang Road, Maogang Town,  
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com



www.sh-bx.com

## Vacuum Oven >>

BZF Series

### Product Feature >

- The shell is made of cold rolled steel plate, the surface is electrostatically sprayed, the inner is made of stainless steel;
- Microcomputer single-chip temperature control system, temperature control, timing, over temperature alarm functions ;
- LCD Display;
- Timing and timing functions ;
- Temperature can be controlled and displayed independently for each heating shelf;
- The door can be adjusted, a silicone rubber door seal, toughened, bulletproof double glass door;
- Power supply recovery function, data will not be lost when power failure and system halted;
- Vacuum pump is optional.



Model	BZF-30	BZF-50	BZF-100	BZF-250	BZF-500
Volume [L]	30	53	100	250	500
Temp. Range [°C]	RT+10~250	RT+10~250	RT+10~250	RT+10~250	RT+10~250
Temp. Fluctuation at 105°C [°C]	±1	±1	±1	±1	±1
Vacuum Degree [Pa]	133	133	133	133	133
Power [W]	450	1200	1400	2800	4000
Power Supply [V, Hz]	220±10, 50/60	220±10, 50/60	220±10, 50/60	220±10, 50/60	380±10, 50
Timer Range [min]	0~9999	0~9999	0~9999	0~9999	0~9999
Internal Dimensions W×D×H [mm]	320×320×300	415×370×345	450×500×450	600×600×700	630×810×845
Housing Dimensions W×D×H [mm]	480×530×615	575×550×675	620×640×780	770×780×1030	800×990×1175
Number of Shelves [pc]	1	2	2	3	4
Net Weight [kg]	45	72	88	175	302
Gross Weight [kg]	55	84	103	194	335
Package Size [mm]	610*580*800	645*690*835	740*770*935	880*865*1243	1090*970*1405